ARMSURE

AM0610 - AM0640

INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 07/2022

Revision: B

EC REP

EC Representative

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Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

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






ArmSure + Shield (AM0610)

	C	Arm Support XR x 2
	D	X-Ray Shield x 1
	B	Arm Strap x 10

Setup

(Instalación / Installation / Configurazione / Einrichtung / Installatie)

X-Ray Shield* (*ArmSure + Shield)

Spare Parts

(Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)

Product Code

(Código de producto / Code produit / Codice prodotto / Produktcode / Productcode)

	ArmSure	AM0640
	ArmSure + Shield	AM0610

Component Codes

(Códigos de componentes / Codes des composants / Codici dei componenti / Komponentencodes / Onderdeelcodes)

	A	Arm Support	M5109
	B	Arm Strap x 10	AM0620
	C	Arm Support XR	M5101
	D	X-Ray Shield	AM0630

Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità tabella immagini / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria della tabella rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

Table Width: / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	45 cm - 70 cm
Table Length: / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	≥200 cm
Table Angulation: / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 15° lateral & longitudinal +/- 15° lateral y longitudinal +/- 15° latéral et longitudinal +/- 15° laterale e longitudinale +/- 15° lateral und longitudinal +/- 15° lateraal en longitudinaal
Table Surface: / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat Plana Plane Piatto Flach Vlak

ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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Information

Intended Purpose

ArmSure: To support the patient's arm.

Arm Strap: To immobilise the patient's arm, protecting the sterile field.

X-Ray Shield: To provide scatter radiation protection for clinician..

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure product is fully inserted under the mattress.
- Do not use with patients over 135Kg.
- Not for paediatric use.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on:

adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No visible damage or sharp edges, e.g. cracks.
- Arm Strap is not damaged, e.g. torn.
- Ensure correct placement of the Arm Support (page 6).
- Ensure correct application of the Arm Straps (page 6-7).

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.
- Any exposed carbon fibre is visible.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

Arm Support/Arm Support XR: Remove Arm Support from beneath the mattress.

Arm Straps: Remove Arm Strap from Arm Support.

X-Ray Shield (when applicable): Remove from Arm Support.

adeptmedicaltraining.com/productcleaning/ArmSure

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

Arm Support, Arm Support XR, Arm Strap and X-Ray Shield

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories
70% Isopropyl Alcohol

Jaychem Industries
2% Chlorhexidine / 70% Alcohol

Betadine
Povidone Iodine (7.5%)

BODE Chemie GmbH
Kohrsolin FF

STERIS
Coverage Spray TB

CaviCide
Metrex Research

Dissolvable:

Du Pont
Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach

PDI Sani-Cloth AF3

Reynard Health Supplies
Surface Disinfectant Wipes

Metrex Research
CaviWipes

Clinell Wipes - Universal (Green)

Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
<5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via adeptmedicaltraining.com/downloads

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

adeptmedical@adept.co.nz

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

X-Ray Shield:

- Ensure no exposed Carbon Fibre is visible.

All components:

- Ensure no visible damage or sharp edges, e.g. cracks.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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