LOWER LEG SUPPORT

AM3000

INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 12/2021

Revision: B

EC REP

EC Representative

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

CONTENTS

Contenido / Contenu / Sommario / Inhalt / Inhoud

- 4** **Unboxing** (Desembalaje / Déballage / Estrazione dalla confezione / Entpacken / Inhoud van de doos)
- 5** **Assembly** (Montaje / Assemblage / Assemblaggio / Montage / Montage)
- 6** **Setup** (Instalación/ Installation / Configurazione / Einrichtung / Installatie)
- 9** **Spare Parts** (Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)
- 10** **Image Table Compatibility** (Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

11 **Essential Information (EN)**

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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


	A	Foot Plate x1
	B	Leg Pad x1
	C	Lower Leg Support x1
	D	Leg Strap x2
	E	Foot Strap x1
	F	Table Strap x2

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




Setup

(Instalación / Installation / Configurazione / Einrichtung / Installatie)

Tightening table straps for device immobilisation.

- / Apriete de las correas de pierna para la inmovilización del dispositivo.
- / Serrage des sangles de table pour l'immobilisation du dispositif.
- / Serraggio delle cinture per il tavolo (Table strap) per l'immobilizzazione del dispositivo.
- / Ziehen Sie die Table Straps fest, um den Leg Support zu fixieren.
- / Tafelriemen (table straps) voor het immobiliseren van het hulpmiddel.

Foot Plate Removal

- / Retirada de la Placa de pie (Foot Plate)
- / Retrait de la plaque de pied (Foot Plate)
- / Rimozione della Foot Plate
- / Entfernen der Foot Plate
- / Verwijdering van voetplaat (Foot Plate)

Spare Parts

(Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)

Product Code

(Código de producto / Code produit / Codice prodotto / Produktcode / Productcode)

	Lower Leg Support	AM3000
---	-------------------	--------

Component Codes

(Códigos de componentes / Codes des composants / Codici dei componenti / Komponentencodes / Onderdeelcodes)

	A	Foot Plate	M6511
	B	Leg Pad	M6512
	C	Support Base	M6516
	D	Leg Strap	M6513
	E	Foot Strap	M6514
	F	Table Strap	M6515

Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilitat Bildgebungstische / Compatibiliteit van MRI-tafel)

Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria del tavolo per imaging rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

Table Width: / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	45cm - 70cm
Table Length: / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	≥200cm
Table Angulation: / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 15° lateral & longitudinal +/- 15° lateral y longitudinal +/- 15° latéral et longitudinal +/- 15° laterale e longitudinale +/- 15° lateral und longitudinal +/- 15° lateraal en longitudinaal
Table Surface: / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat Plana Plane Piatto Flach Vlak

ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Wichtige Informationen / Essentiële informatie

11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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Information

Intended Purpose

To support and immobilise a supine patient's leg during endovascular surgery.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not over tighten Leg & Foot Straps.
- Do not use with patients over 135Kg.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on

adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks.
- **Straps:** ensure there is no damage, e.g. tearing.
- **Table Strap:** ensure buckle can fully engage and disengage smoothly.
- **Table Strap:** ensure Table Strap is secure and correctly placed (page 6).
- **Foot/Leg Strap:** ensure straps are secure and correctly placed (page 7 and 8).
- **Foot Plate:** Ensure the Foot Plate is correctly placed in the Foot Plate Slot and fully engaged.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.
- Any exposed carbon fibre is visible.
- The Foot Plate does not securely fit into the Lower Leg Support Base.
- Any tearing is observed on the Leg Pad or any of the Straps.
- The hook & loop on the Leg/Foot Strap does not securely engage.
- The buckle on the Table Strap does not securely engage.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

Lower Leg Support Base: Remove the base from the mattress.

Foot Plate: Detach the Foot Plate from the Base.

Leg Pad: Remove the padding from the Base.

Leg Strap: Unfasten the strap i.e. ensure the "hook" and "loop" sides are not connected.

Foot Strap: Unfasten the strap i.e. ensure the "hook" and "loop" sides are not connected.

Table Strap: Unfasten the buckle on the strap.

adeptmedicaltraining.com/productcleaning/Lower-Leg-Support

Cleaning: Manual

Lower Leg Support Base/Foot Plate

- If visibly soiled rinse with water.
- These components can be fully submerged in water.
- A soft bristled brush can be used if required.
- Ensure soiling is fully removed.
- Ensure all components are completely dry prior to disinfection.

Leg Pad/Leg Strap/Foot Strap/Table Strap

- If visibly soiled rinse with water.
- These components should not be fully immersed in water.
- A soft bristle brush can be used if required.
- Ensure soiling is fully removed.
- Ensure all components are fully dry prior to disinfection.

Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

Lower Leg Support:

With a cloth dampened by an approved disinfectant, wipe all surfaces, including inside the Foot Plate slot.

Foot Plate/Leg Pad/Leg Strap/Foot Strap and Table Strap:

With a cloth dampened by an approved disinfectant, wipe all surfaces.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

- ORION Laboratories
- 70% Isopropyl Alcohol

- Jaychem Industries
- 2% Chlorhexidine / 70% Alcohol

- Betadine
- Povidone Iodine (7.5%)

- BODE Chemie GmbH
- Kohrsolin FF

- STERIS
- Coverage Spray TB

- CaviCide
- Metrex Research

Dissolvable:

- Du Pont
- Rely+On Virkon Tablets

Surface Wipes:

- PDI Sani-Cloth Bleach

- PDI Sani-Cloth AF3

- Reynard Health Supplies
- Surface Disinfectant Wipes

- Metrex Research
- CaviWipes

- Clinell Wipes - Universal (Green)

- Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

- ≤10% Malic Acid CAS 6915-15-7
- ≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

- ≤5% 2-Butoxyethanol CAS 111-76-2
- ≤10% Butyldiglycol CAS 112-34-5
- ≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
- ≤10% Tridecanol CAS 69011-36-5
- <10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
- ≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

- ≤0.1% PHMB CAS 27083-27-8
- ≤55% Potassium Peroxymonosulfate CAS 70693-62-8
- ≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

- ≤5% Benzalkonium Chloride CAS 68424-85-1
- ≤0.28% Benzethonium Chloride CAS 121-54-0
- ≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
- ≤2% Chlorhexidine CAS 55-56-1
- ≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
- ≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
- <5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

Others:

- ≤10% Glutaral CAS 111-30-8
- ≤0.63% Sodium Hypochlorite CAS 7681-52-9
- ≤7.5% Povidone Iodine CAS 25655-41-8
- 0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
- ≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via adeptmedicaltraining.com/downloads

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

adeptmedical@adept.co.nz

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

Lower Leg Support/Foot Plate:

- No exposed Carbon Fibre is visible.
- No visible damage or sharp edges, e.g. cracks.
- Ensure that the Foot Plate securely fits in to Lower Leg Support.

Leg Pad:

- No visible damage or tearing.

Foot Strap/Leg Strap:

- No visible damage or tearing.
- The hook & loop should attach securely.

Table Strap:

- No visible damage or tearing.
- The buckle can fully engage and disengage smoothly.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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