MR Safe
**OVERHEAD
ARM SUPPORT**
AM4000

INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 07/2022

Revision: B

EC Representative

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



	A	Armrest Wing x1 MR Safe
	B	Arm Strap x 4
	C	Locking Leg Assembly x1 MR Safe
	D	Daggerboard x1 MR Safe

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








	A	Armrest Wing MR Safe	M5504
	B	Arm Strap x 10	AM0620
	C	Locking Leg Assembly MR Safe	M5505
	D	Daggerboard MR Safe	M5506

Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità tabella immagini / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Assicurarsi che la geometria della tabella rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

Table Width: / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	45cm - 70cm
Table Length: / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	≥200cm
Table Angulation: / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 15° lateral & longitudinal +/- 15° lateral y longitudinal +/- 15° latéral et longitudinal +/- 15° laterale e longitudinale +/- 15° lateral und longitudinal +/- 15° lateraal & longitudinaal
Table Surface: / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat / Curved Plana / Curva Plane / incurvée Piatto / Curvo Flach / Gebogen Vlak / gebogen
CT/MRI Bore: / Diámetro del túnel de TC/RM / Alésage de l'appareil de TDM/IRM / Diametro interno CT/MRI / Tunneldurchmesser CT/MRT / CT/MRI-opening	≥70cm

ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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Information

Intended Purpose

To support the patient's arms overhead with controlled shoulder flexion for imaging. Intended for bore's 70cm and upward.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.
- Not intended for bore size less than 70cm.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure MRI safety checks are conducted.
- Ensure patient and device will fit clear of the bore.
- Do not use with patients over 135Kg.
- Not for paediatric use.
- Ensure general patient insulation protocol is followed for MRI.
- For CT and MRI, ensure there is adequate patient clearance prior to use. (refer to page 8)

Incident Reporting Guidelines

For product complaints and incidents, please complete form on:

adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.


Product Checks

- No visible damage or sharp edges, e.g. cracks.
- **Arm Strap:** not damaged, e.g. tearing.
- **Locking Leg Assembly:** ensure that the Locking Leg Assembly is not free to move in the locked position.
- **Locking Leg Assembly:** Ensure it can move up and down when pushed in the unlocked position.
- Ensure the Locking Leg Assembly is fully inserted into the Daggerboard and the release clip is fully engaged (refer to page 5).
- Ensure the Armrest Wing is fully inserted onto the Locking Leg Assembly (refer to page 5).


MRI Safety Checks

Ensure each component below has the 'MR Safe' product sticker and engraving in the below locations:


Armrest Wing:


Locking Leg Assembly:

Daggerboard:

Safety Clearance refer to page 8

Patient Management

Ensure insulation protocol is followed according to individual facility guidelines. This includes ensuring sufficient insulation is placed to eliminate skin to skin conductive loops and that the patient's skin does not contact the bore of the magnet, when using the device.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

- Remove Overhead Arm Support from beneath the mattress.
- Remove Arm Straps from Armrest Wing.
- Disassemble Overhead Arm Support into Armrest Wing, Locking Leg Assembly and Daggerboard.
- Ensure the Locking Leg Assembly is fully extended.

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

Armrest Wing

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including inside the Arm Strap slots.

Daggerboard

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including inside of the Locking Leg Assembly slot.

Locking Leg Assembly

- With a cloth dampened by an approved disinfectant, wipe all surfaces, especially the blue lever.

Arm Strap

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

adeptmedicaltraining.com/productcleaning/Overhead-Arm-Support

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

- ORION Laboratories
- 70% Isopropyl Alcohol

- Jaychem Industries
- 2% Chlorhexidine / 70% Alcohol

- Betadine
- Povidone Iodine (7.5%)

- BODE Chemie GmbH
- Kohrsolin FF

- STERIS
- Coverage Spray TB

- CaviCide
- Metrex Research

Dissolvable:

- Du Pont
- Rely+On Virkon Tablets

Surface Wipes:

- PDI Sani-Cloth Bleach

- PDI Sani-Cloth AF3

- Reynard Health Supplies
- Surface Disinfectant Wipes

- Metrex Research
- CaviWipes

- Clinell Wipes - Universal (Green)

- Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

- ≤10% Malic Acid CAS 6915-15-7
- ≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

- ≤5% 2-Butoxyethanol CAS 111-76-2
- ≤10% Butyldiglycol CAS 112-34-5
- ≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
- ≤10% Tridecanol CAS 69011-36-5
- <10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
- ≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

- ≤0.1% PHMB CAS 27083-27-8
- ≤55% Potassium Peroxymonosulfate CAS 70693-62-8
- ≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

- ≤5% Benzalkonium Chloride CAS 68424-85-1
- ≤0.28% Benzethonium Chloride CAS 121-54-0
- ≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
- ≤2% Chlorhexidine CAS 55-56-1
- ≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
- ≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
- <5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

Others:

- ≤10% Glutaral CAS 111-30-8
- ≤0.63% Sodium Hypochlorite CAS 7681-52-9
- ≤7.5% Povidone Iodine CAS 25655-41-8
- 0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
- ≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via adeptmedicaltraining.com/downloads

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

adeptmedical@adept.co.nz

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

For all components ensure:

- No visible damage or sharp edges, e.g. cracks.

Locking Leg Assembly:

- Ensure blue lever functions smoothly.
- Ensure when lever is in the unlocked position the Locking Leg Assembly can extend and contract smoothly when force is applied.
- Ensure when lever is in the locked position the Locking Leg Assembly cannot contract or extend when force is applied.
- Ensure smooth engagement with daggerboard during assembly.
- Ensure tight fit with Armrest Wing during assembly.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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