

# Hega

## *Dina & Nova Rinsing Guns, and Attachments*



### **Medical Resources Group**

1920 South Main Street

Morton, IL 61550

+1 309.263.6901 | [sales@mrg-ltd.com](mailto:sales@mrg-ltd.com) | [www.mrg-ltd.com](http://www.mrg-ltd.com)



# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5




309.263.6901

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)

Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
<b>Dina - Rinsing Gun</b>			
<b>Dina Rinsing Gun without hose and attachments</b>	<b>Connection</b>	<b>Item No.</b>	
Dina Rinsing Gun with Spout		6000	
<b>Dina Rinsing Gun with open end of hose and a mountable connection set</b>	<b>Connection</b>	<b>Item No.</b>	
1,5m PVC-hose, fabric reinforcement	1/2"	6015	
	3/4"	6025	
3,0m PVC-hose, fabric reinforcement	1/2"	6035	
	3/4"	6045	
1,5m stainless steel hose	1/2"	6055	
	3/4"	6065	
3,0m stainless steel hose	1/2"	6075	
	3/4"	6085	
<b>Dina Cleaning-Kit with Professional Set including 8 attachments and an attachment holder</b>	<b>Connection</b>	<b>Item No.</b>	
1,5m PVC-hose, fabric reinforcement	1/2"	6010	
	3/4"	6020	
3,0m PVC-hose, fabric reinforcement	1/2"	6030	
	3/4"	6040	
1,5m stainless steel hose	1/2"	6050	
	3/4"	6060	
3,0m stainless steel hose	1/2"	6070	
	3/4"	6080	

\*Connection Set: 3/4" is a standard garden hose connection



# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)




Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
---------------------	---------	----------	-----------------

## Nova - Rinsing Gun

Nova Rinsing Gun without hose and attachments	Connection	Item No.	
Nova Rinsing Gun with Spout		7000	
Nova Rinsing Gun with open end of hose and a mountable connection set	Connection	Item No.	
1,5m PVC-hose, fabric reinforcement	1/2"	7015	
	3/4"	7025	
3,0m PVC-hose, fabric reinforcement	1/2"	7035	
	3/4"	7045	
1,5m stainless steel hose	1/2"	7055	
	3/4"	7065	
3,0m stainless steel hose	1/2"	7075	
	3/4"	7085	
Nova Cleaning-Kit with Professional Set including 8 attachments and an attachment holder	Connection	Item No.	
1,5m PVC-hose, fabric reinforcement	1/2"	7010	
	3/4"	7020	
3,0m PVC-hose, fabric reinforcement	1/2"	7030	
	3/4"	7040	
1,5m stainless steel hose	1/2"	7050	
	3/4"	7060	
3,0m stainless steel hose	1/2"	7070	
	3/4"	7080	

\*Connection Set: 3/4" is a standard garden hose connection



# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5











309.263.6901

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)

Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
<b>Equipment for the Dina and Nova Rinsing Guns</b>			
Pointy attachment, for the cleaning of valves		4230	
Drainage attachment, for the cleaning of drainage hose		4240	
Spray jet attachment, for the cleaning of surfaces and hollows		4250	
Record attachment with record cone, for the cleaning of syringes and cannulas		4270	
Valve attachment, for the cleaning of valves		5010	
Luer attachment with luer connection, for the cleaning of syringes and cannulas		4260	
Bottle cleaner attachment, for the cleaning of bottles and Erlenmeyer flasks		5020	
Water jet pump attachment, for the suction of angles and hollows		5030	
<b>Professional Set.</b> Including 8 attachments and 1 attachment holder, stainless steel 2mm, 280 x 40mm, self-adhesive		5220	
Attachment holder for 8 attachments, stainless steel 2mm, 280 x 40mm, self-adhesive		5040	

\*Connection Set: 3/4" is a standard garden hose connection