





# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)

Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
---------------------	---------	----------	-----------------

Hose Connection with thread, open end of hose			
crimped with spout (5150+3010) for Rinsing Gun, o-ring 13x2,5 (4551), cone (4390), connection set	Connection	Item No.	
1,5m PVC-hose, fabric reinforcement	1/2"	S600	
	3/4"	S602	
3,0m PVC-hose, fabric reinforcement	1/2"	S610	
	3/4"	S612	
1,5m stainless steel hose	1/2"	S620	
	3/4"	S622	
3,0m stainless steel hose	1/2"	S630	
	3/4"	S632	
Quick-Coupling Spiral Hose		Item No.	
5m Spiral hose for compressed air, black PU, with quick release coupling and plug DN 7.2, nickel-plates steel diameter of hose 6mm/8mm		S150	
Service Kit			
Service kit for Dina-Rinsing Gun. Content: sealing cone 1008, dosing pin 2050, dosing wheel 2060, lever, spiral pin, compression spring 1010, closure lid 2012, o-ring 13x2,5 mm 4551, o-ring 3x1mm 4550		S014	

\*Connection Set: 3/4" is a standard garden hose connection



# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

[sales@mr-g-ltd.com](mailto:sales@mr-g-ltd.com)

Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
<b>Equipment for the Dina and Nova Rinsing Guns</b>			
Connection set 1/2". Including chrome plated brass spout 1/2", coupling nut 1/2", 2 x jurid clamp, silicone gasket	1/2"	5110	
Connection set 3/4". Including chrome plated brass spout 3/4", coupling nut 3/4", 2 x jurid clamp, silicone gasket	3/4"	5120	
Spout for Dina and Nova Rinsing Guns with M16, external thread for cone (4390) and O-Ring 13x2,5mm		5150	
Cone for Dina and Nova Rinsing Guns with M16 internal thread, incl. O-ring 13x2,5mm, fits in table implementation (5070) and wall holders (3000)/(3006)		4390	
Outlet connector for medical air, including 2 jurid clamps for installation	DIN 13260	8140	
Nozzle, for Dina Rinsing Guns, M12x1,25 with O-Ring 13x2,5 and 3x1mm, stainless steel	M12X1,25	5130	
Nozzle, for Nova Rinsing Guns with O-Ring 3x1mm, chromed brass	UNS 1/2" - 27	5140	
Ferrule stainless steel, for PVC or stainless steel hose. Price per piece		3010	
PVC hose, fabric reinforcement. Price per meter		3025	
Stainless steel hose with stainless steel sheath tissue. Price per meter		3030	
Reducer, reduction of 3/4" male to 1/4" female, chromed brass		5160	
Jurid clamp, stainless steel. Price per piece		4520	

\*Connection Set: 3/4" is a standard garden hose connection



# Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)

Effective 1/1/26

[www.mrg-ltd.com](http://www.mrg-ltd.com)

[www.mrg-ltd.com/hega-medical.html](http://www.mrg-ltd.com/hega-medical.html)

Product Description	Details	Item No.	Product picture
<b>Equipment for Dina and Nova Rinsing Guns</b>			
Wall holder for all HEGA Rinsing and Washing Guns, for wall installation. Material: Nylon		3000	
Wall holder for all HEGA Rinsing and Washing Guns, for wall installation. High shine polished stainless steel		3006	
Stainless steel hook with 2 Phillips screws		3020	
Table implementation for all HEGA Rinsing and Washing Guns. With lock nut SW42 and silicone gasket. For table mounting, hole diameter 34 to 40mm, table thickness up to 50mm		5070	
Table implementation, stainless steel, for modified cleaning guns. With lock nut SW42 and silicone gasket. For table mounting, hole diameter 34 to 40mm, table thickness up to 50mm		4482	
Gasket 1/2" silicone washer ID=8mm, OD=19mm. Price per piece	1/2"	4530	
Gasket 3/4" silicone washer ID=8mm, OD=25mm. Price per piece	3/4"	4540	
O-Ring 13x2,5mm, silicone. Price pro Stück (Pack of 3)	Modell: Dina	4551-3	
O-Ring for Dina- Nova- und Hela- Rinsing Guns, silicone 3x1mm (Pack of 3)		4550-3	

\*Connection Set: 3/4" is a standard garden hose connection