

Hega

Spare Parts for Hela Rinsing Guns



Medical Resources Group

1920 South Main Street

Morton, IL 61550

+1 309.263.6901 | sales@mrg-ltd.com | www.mrg-ltd.com



Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

sales@mrg-ltd.com

Effective 1/1/26

www.mrg-ltd.com

www.mrg-ltd.com/hega-medical.html

Product Description	Details	Item No.	Product picture
Spare parts and equipment for the Hela Stainless Steel Cleaning Gun			
Dosing lever for the Hela Cleaning Gun		9032	
Silicone handle for the Hela Cleaning Gun		9042	
Spout for the Hela Cleaning Gun		1030	
Holding cone for the Hela Cleaning Gun, incl. O-ring 13x2mm, fits in table implementation (5070) and wall holders (3000)/ (3006).		1035	
Universal nozzle for the Hela Cleaning Gun incl. O-ring 13x2mm to fit the attachments of the Professional Set (5220) and Brush Set (5210)		1100	
Wall holder for all HEGA Rinsing and Washing Guns, for wall installation. Material: Nylon		3000	
Wall holder for all HEGA Rinsing and Washing Guns, for wall installation. High shine polished stainless steel.		3006	
Table implementation for all HEGA Rinsing and Washing Guns. With lock nut SW42 and silicone gasket. For table mounting, hole diameter 34 to 40mm, table thickness up to 50mm		5070	
Table implementation, stainless steel, for modified cleaning guns. With lock nut SW42 and silicone gasket. For table mounting, hole diameter 34 to 40mm, table thickness up to 50mm		4482	
O-Ring for Dina- Nova- und Hela- Rinsing Gun, silicone 3x1mm (Pack of 3)		4550-3	
O-Ring for the Hela Cleaning Gun, silicone 13x2mm (Pack of 3)		4560-3	

*Connection Set: 3/4" is a standard garden hose connection



Medical Resources Group

1920 South Main Street Morton, IL 61550

Price List V9.5

309.263.6901

sales@mrg-ltd.com

Effective 1/1/26

www.mrg-ltd.com

www.mrg-ltd.com/hega-medical.html

Product Description	Details	Item No.	Product picture
Hela Hose Connection with thread, open end of hose			
crimped with spout (1030+3010) for Hela Rinsing Gun, o-ring 13x2 (4560), cone (1035), connection set	Connection	Item No.	
1,5m PVC-hose, fabric reinforcement	1/2"	S900	
	3/4"	S902	
3,0m PVC-hose, fabric reinforcement	1/2"	S910	
	3/4"	S912	
1,5m stainless steel hose	1/2"	S920	
	3/4"	S922	
3,0m stainless steel hose	1/2"	S930	
	3/4"	S932	
Quick coupling spiral hose		Item No.	
5m Spiral hose for compressed air, black PU, with quick release coupling and plug DN 7.2, nickel-plated steel diameter of hose 6mm/8mm		S150	
Service Kit for Hela			
Service Kit for Hela Cleaning Gun. Content: sealing cap 1025, compression spring 1011, conical nipple 1008, dosing pin 1020, trigger axle 1015, O-ring 13x2mm 4560		S016	

*Connection Set: 3/4" is a standard garden hose connection