


# CaraClip™

wire protector

**Beekley Medical®**  
when your diagnosis  
must be right®


**Reduce risk of wire contamination, free up hands,  
and organize wires in the sterile field.**




CaraClip™ wire protector was designed by an interventional radiologist to firmly adhere to the drape and contain the guidewire and other accessories, unlike unreliable methods like towels and gauze.

- Helps prevent contamination by keeping guidewires and other accessories in the sterile field.
- Helps curve, redirect, and allows for the advancement of guidewires, catheters, tubing/sheaths, and other accessories.
- Remains in place for at least 8 hours on dry drape and withstands fluids such as: blood, saline, iodine, contrast, and sterile water.
- Repositionable on dry drape and remains in the field for at least 8 hours.
- Helps manage guidewires and other accessories on the patient drape and procedure table, saving time and aiding in organization.

### Back table organization:



Guidewire coiled and secured by CaraClip with Telfa™



Guidewire coiled and secured by CaraClip



REF 431

QTY 15 pouches of 2 units

STERILE R

### Distributed By:



### Medical Resources Group


1920 South Main Street

Morton, IL 61550

+1-309-263-6901




[www.mrg-ltd.com](http://www.mrg-ltd.com)

[sales@mrg-ltd.com](mailto:sales@mrg-ltd.com)

 Manufactured by Beekley Corporation  
One Prestige Lane, Bristol, CT 06010-7454 USA  
1.800.233.5539 or +1.860.583.4700 [beekley.com](http://beekley.com)

CARA CLIP is an unregistered trademark of Beekley Corporation.  
All other trademarks are the property of their respective owners.  
® – Reg. U.S. Pat. & Tm. Off.  
© 2023 Beekley Corporation. All rights reserved.  
Pat. – [beekley.com/patents](http://beekley.com/patents)

### WARNINGS AND PRECAUTIONS:

- Do not use if product sterile packaging is compromised. 
- Single patient use only; Do not reuse. 
- Do not resterilize. 
- Does not contain natural rubber latex.



REV: CC\_INF\_1023