

## Stainless Steel Closed Case Carts



### Standard features:

- Welded 304 stainless steel heavy duty construction
- High temperature heavy duty bumper-will not shrink or hold water
- Sealed insulated heavy duty double walled doors with transport lock out latch
- Ergonomic push handles
- Sealed heavy duty no maintenance casters
- Positive lock door latches
- Stainless steel adjustable wire shelf(s) with smooth fronts

### Models:

---

999EHW1S1DT0	1 Door Short	33" W x 30.3" D x 39.4" H
999EHW1S2DT0	2 Door Short	48.8" W x 30.3" D x 39.4" H
999EHW2S1DT1	1 Door Tall	30.9" W x 29.5" D x 57.8" H
999EHW2S2DT0	2 Door Tall	46.6" W x 29.5" D x 57.8" H

### **STEELCO USA Inc.**

7354 Central Industrial Drive West Palm Beach FL 33404 USA

Ph: 561 791 8313 Fax 561 660 5382 www.steelcospa.com e-mail: contact@steelco-usa.com

Optional Items:

- Solid stainless steel adjustable shelves
- Paper document clamp
- Pull out sliding shelves
- Glass doors
- Laser barcode ID plate
- 3 or 5 Color Wheel
- Clean/Hazardous Sign

