

ID300 Drying Cabinet for Instruments



Models:

ID300-1	single door	27.95" W x 30.71" D x 74.8" H
ID300-2	pass-through	27.95" W x 32.48" D x 74.8" H

Standard Features:

- HEPA filtration of heated drying air
- UV-light disinfection of air flow circuit
- LED Display
- Supplied with 8 Shelves
- "Air curtain" assuring overpressure in cabinet with opened door
- Closed stainless steel cabinet with tempered glass lockable swing door(s)
- Temperature, air flow, and door lock alarms
- Completely Self Contained

Options:

- Humidity Sensor
- Seismic Anchoring Kit
- 8 Luer Lock Connections for cannulated instruments

Available Configurations:

- 120 V, 60Hz single phase

STEELCO USA Inc.

7354 Central Industrial Drive West Palm Beach FL 33404 USA

Ph: 561 791 8313 Fax 561 660 5382 www.steelcospa.com e-mail: contact@steelco-usa.com