

# Hygienic Covers for Prone Support Face Pad



Align the Velcro strip in the center of the  
Prone Support Face Pad



Attach the top and bottom Velcro strips centered  
on the back side of the Face Pad



Attached Prone Support Face Pad to the  
Prone Support

Product Number: MRG0700FC

Pack / 50

Available From:  
**Medical Resources Group**  
1920 South Main Street Morton, IL 61550  
800-674-6621 | [sales@mrg-ltd.com](mailto:sales@mrg-ltd.com) | [www.mrg-ltd.com](http://www.mrg-ltd.com)